
Kraftgürtel

aus recyceltem Kunststoff

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Kurt-Georg-Kiesinger-Allee 3 - 53175 Bonn

<https://www2.bfarm.de/medprod/mpsv/>

GB | Support belt, made from recycled plastic

1 | Delivery

1 x Support belt

Read the instructions for use completely and carefully!

2 | Intended purpose

By external pressure (compression), the abdominal and back muscles should be activated and supported. Supports the spine and promotes an upright posture.

- Disc problems
- Pain of sciatica
- Lumbago
- Pain from overstrain

For private use

3 | Contraindications / Not suitable ⚠

In the case of the following medical conditions, it is only recommended to wear this article after consultation with your doctor:

- Skin diseases/injuries in the affected part of the body, especially in inflammatory symptoms such as swelling, redness and overheating
- Allergic reactions to the material
- Sensory and movement disorders, e.g. in the case of diabetes
- Open wounds or internal injuries
- For small children or persons who can not be noticeable when the product sits uncomfortably

4 | Notes on safe use of the product ⚠

Side effects affecting the whole organism are not yet known. Proper application is required. In case of improper use, product liability is excluded.

- If the bandage is too tight or uncomfortable, remove it immediately!
- Used by the times by private users.

- If your symptoms worsen or you notice unusual changes, consult your doctor immediately.

Do not use:

While sleeping

5 | Application

Put on the support belt: Fig.1: Hold the support belt with the hook and loop closure (inside) in the right hand and guide it around the back. Pull the two straps sideways away from the body. Fig.2: First fold in the left side and then the right one. The first closure allows a simplified system. Fig.3: Then close the hook and loop fastener at any point on the bandage.

6 | Cleaning

Close hook and loop fasteners during washing to avoid damaging other items of laundry. Cleaning is done by household laundry, according to care instructions.

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